



CL016 English

### CL016 Mini Vortex

Suitable for mixing the contents of specimen tubes in the laboratory.  
The instrument operates automatically upon pressing on the head of the vortex (press start function).

**Choose **CYAN CL016** and get the following advantages:**

- Automatic press start function
- Variable speed control
- Compact, robust, and durable
- Vibration free
- Fast
- Quiet
- High-quality motor for perfect performance

**Technical data:**

- Silicon cup head
- Orbit of 4,5 mm
- Max. speed: 3200 rpm
- Max. volume: 50 ml
- Max. weight: 60 g
- Head cup dimensions: 3,9 cm (Ø)
- Overall dimensions: 5,5 x 11,3 x 8,1 cm (W x D x H)
- Gross weight: 0,90 kg
- Power supply: 110/220 V - 50/60 Hz

**Order code:**

- CL016: **Mini Vortex** (110/220 V - 50/60 Hz)

# Mini Vortex



Contact: e-mail: [cypress@diagnostics.be](mailto:cypress@diagnostics.be)  
 tel: + 32 (0)15 67 67 68  
 Nijverheidsstraat 8  
 2235 Hulshout  
 Belgium

URL: [www.diagnostics.be](http://www.diagnostics.be)

Your Distributor

Mini Vortex EN 2019-04