



CL010 English

CL010 Water bath

Choose **CYAN CL010** and get the following advantages:

- Low electricity consumption
- Very stable temperature
- Bath with rounded corners for easier cleaning

Technical Data:

- Unstirred (thermal diffusion)
- Working range: ambient temperature up to 100 °C
- Capacity: 10 liters
- PID temperature control
- Digital temperature setting
- Digital temperature readout
- Stability: +/- 0,2 °C
- Inner tank dimensions: 24 x 30 x 15 cm (W x D x H)
- Overall dimensions: 29 x 40 x 24 cm (W x D x H)
- Gross weight: 8,5 kg
- Overheating protection
- Power: 650 W
- Power supply: 110 V – 60 Hz (CL010-110)
or 220 V – 50/60 Hz (CL010-220)

Accessories included:

- Flat cover for water bath

Order code:

- CL010-110: [Water bath](#) (110 V – 60 Hz)
- CL010-220: [Water bath](#) (220 V – 50/60 Hz)

Auxiliary products:

- CL010-R Test tube rack 50 places

Water Bath



Contact: e-mail: cypress@diagnostics.be
tel: + 32 (0)15 67 67 68
Nijverheidsstraat 8
2235 Hulshout
Belgium
URL: www.diagnostics.be

Your Distributor

Waterbath EN 2019-04



ISO 13485-2016

www.diagnostics.be • Belgium • Tel: ++ 32 15 67 67 68 • e-mail: cypress@diagnostics.be



ISO 13485-2016

www.diagnostics.be • Belgium • Tel: ++ 32 15 67 67 68 • e-mail: cypress@diagnostics.be