

Control blood for CYANHemato/CYANHemato 380							
		Control LOW		Control NORMAL		Control HIGH	
		Lot: 61520422		Lot: 61520423		Lot: 61520424	
		Exp: 2026-12-07		Exp: 2026-12-07		Exp: 2026-12-07	
		Open vial stability: 30 days		Open vial stability: 30 days		Open vial stability: 30 days	
Parameter	Unit	Target	Range (+/-)	Target	Range (+/-)	Target	Range (+/-)
WBC	10 ⁹ /l	2,50	0,40	7,30	0,60	20,30	2,50
RBC	10 ¹² /l	2,31	0,15	4,21	0,20	5,20	0,25
HGB	g/dl	5,4	0,5	11,2	0,6	15,6	0,7
[HGB]	g/l	54	5	112	6	156	7
HGB	mmol/l	3,4	0,3	7,0	0,4	9,7	0,4
HCT	%	17,10	2,00	34,00	3,00	44,40	4,00
[HCT]	l/l	0,171	0,020	0,340	0,030	0,444	0,040
MCV	fl	74,8	6,0	80,2	6,0	86,1	6,0
MCH	pg	22,9	2,0	26,4	2,0	29,9	2,0
MCH	fmol	1,45	0,12	1,65	0,12	1,86	0,12
MCHC	g/dl	32,0	3,0	33,4	3,0	35,1	3,0
[MCHC]	g/l	320	30	334	30	351	30
MCHC	mmol/l	19,6	0,2	20,4	0,2	21,8	0,2
PLT	10 ⁹ /l	75	15	226	30	636	60
PCT	%	0,05	0,03	0,17	0,06	0,47	0,20
MPV	fl	7,6	1,5	7,5	1,5	7,9	1,5
PDWc	%	34,4	5,5	32,2	3,7	33,0	2,7
RDWc	%	18,5	5,0	18,3	5,0	17,0	5,0
LYM	10 ⁹ /l	1,50	0,50	3,00	0,50	4,00	1,50
MID	10 ⁹ /l	0,20	0,20	0,60	0,30	1,90	1,00
GRA	10 ⁹ /l	0,80	0,80	3,70	1,00	14,40	3,00
LYM	%	60,1	8,0	41,1	6,0	19,7	6,0
MID	%	8,0	6,0	8,2	5,0	9,3	5,0
GRA	%	31,9	8,0	50,7	7,0	71,0	8,0

Product code	Item
CY006-C01A	Control Package Normal (1 x 2.5 ml)
CY006-C02A	Control Package Low/Normal/High (3 x 2.5 ml)

PAY ATTENTION:

before closing the control blood vial, please remove the blood from the cap using lint-free tissue



Use the **lot-specific .ini files** for quick and easy updating of the control values on your CYANHemato/CYANHemato 380.
You can find them at:
<https://diagnostics.be/update/CY006-C0>

Sang de contrôle pour CYANHemato/CYANHemato 380							
		Contrôle BAS		Contrôle NORMAL		Contrôle HAUT	
		Lot: 61520422		Lot: 61520423		Lot: 61520424	
		Exp: 2026-12-07		Exp: 2026-12-07		Exp: 2026-12-07	
		Stabilité après ouverture: 30 jours		Stabilité après ouverture: 30 jours		Stabilité après ouverture: 30 jours	
Paramètre	Unité	Cible	Gamme (+/-)	Cible	Gamme (+/-)	Cible	Gamme (+/-)
GB	10 ⁹ /l	2,50	0,40	7,30	0,60	20,30	2,50
GR	10 ¹² /l	2,31	0,15	4,21	0,20	5,20	0,25
HB	g/dl	5,4	0,5	11,2	0,6	15,6	0,7
[HB]	g/l	54	5	112	6	156	7
HB	mmol/l	3,4	0,3	7,0	0,4	9,7	0,4
HT	%	17,10	2,00	34,00	3,00	44,40	4,00
[HT]	l/l	0,171	0,020	0,340	0,030	0,444	0,040
VGM	fl	74,8	6,0	80,2	6,0	86,1	6,0
TCMH	pg	22,9	2,0	26,4	2,0	29,9	2,0
TCMH	fmol	1,45	0,12	1,65	0,12	1,86	0,12
CCMH	g/dl	32,0	3,0	33,4	3,0	35,1	3,0
[CCMH]	g/l	320	30	334	30	351	30
CCMH	mmol/l	19,6	0,2	20,4	0,2	21,8	0,2
PLT	10 ⁹ /l	75	15	226	30	636	60
THT	%	0,05	0,03	0,17	0,06	0,47	0,20
VPM	fl	7,6	1,5	7,5	1,5	7,9	1,5
IDPc	%	34,4	5,5	32,2	3,7	33,0	2,7
IDRc	%	18,5	5,0	18,3	5,0	17,0	5,0
LYM	10 ⁹ /l	1,50	0,50	3,00	0,50	4,00	1,50
MID	10 ⁹ /l	0,20	0,20	0,60	0,30	1,90	1,00
GRA	10 ⁹ /l	0,80	0,80	3,70	1,00	14,40	3,00
LYM	%	60,1	8,0	41,1	6,0	19,7	6,0
MID	%	8,0	6,0	8,2	5,0	9,3	5,0
GRA	%	31,9	8,0	50,7	7,0	71,0	8,0

Code de produit	Article
CY006-C01A	Coffret de Contrôles Normaux (1 x 2.5 ml)
CY006-C02A	Coffret de Contrôles Bas/Normaux/Hauts (3 x 2.5 ml)

ATTENTION:

avant de fermer le tube du sang de contrôle, enlever le sang du bouchon utilisant du tissu non pelucheux



Utilisez **les fichiers .ini par lot** pour une mise à jour rapide et facile des valeurs de contrôle sur votre CYANHemato (380).
Vous pouvez les trouver à l'adresse suivante :
<https://diagnostics.be/fr/update/CY006-C0>